

Takes the pressure



VERTEBRIS stenosis

Full-endoscopic, interlaminar decompression
in case of lumbar spinal canal stenosis

VERTEBRIS *stenosis*

Full-endoscopic Spine Instrumentation



Contents

Full-endoscopic, interlaminar decompression	
■ Introduction	04
■ Positioning	06
■ Determination of interlaminar access	06
■ Performance of interlaminar access	07
■ Ipsilateral, decompression on one side	08
■ Contralateral decompression in over-the-top technique	10
VERTEBRIS stenosis Instrumentation	
■ Endoscope and accessories VERTEBRIS stenosis	12
■ Access and working instruments VERTEBRIS stenosis	13
■ PowerDrive ART1 – Universal motor system	14
■ COMBIDRIVE EN – High-speed motor system	15
■ Surgitron – Radio frequency unit	16
■ Consumables and replacement parts	17
Literature	18

VERTEBRIS stenosis



Introduction

Degenerative stenosis of the lumbar spinal canal with compression of neural elements arise as a result of bony, disk, capsular or ligament structures. Depending on localization and spread, they can lead to classic symptoms in the lower extremities. Pain in the back tends to be attributed to secondary degenerative phenomena, e.g. segmental instability or deformity. There is no clear correlation between the extent of the stenosis shown by imaging and the clinical symptoms. Apart from spinal disk herniations, lateral and central spinal canal stenoses form the most frequent causes.

A surgical procedure may be necessary after conservative measures have been exhausted or neurological deficits occur. When this is the case, the pathology and symptoms must be taken into account and decompression operations, fusions, or a combination of both procedures must be considered. Today, EBM criteria appear to provide certainty that decompression procedures can improve radicular symptoms and neurogenic claudication. The extent of decompression required from a technical perspective and the conditions under which an additional fusion is necessary have not been definitively described.

Conventional decompression operations on the lumbar spine demonstrate good results. However, consequences and problems associated with these operations are known. Attempts were therefore made right from the start of spine surgery to modify existing operating procedures. Up to the present day, the primary focus continues to be on reducing the invasiveness of surgery and improving the intraoperative view.

Minimally invasive techniques can reduce the trauma and consequences due to the operation. At the same time, visualization and illumination during the operation can be optimized. Appropriate instrument sets for decompression of lumbar spinal canal stenosis were developed on the basis of experiences derived from full-endoscopic operations on spinal disk herniations of the cervical and lumbar spine, offering the possibility of endoscopic bone resection. Since a more extensive bone or ligament resection is frequently necessary here, a large endoscope with a correspondingly large intraendoscopic working channel and larger instruments were necessary. Full-endoscopic, interlaminar access is used routinely, while the transforaminal/extraforaminal access is reserved for specific individual cases.



Lateral and central spinal canal stenosis of the lumbar spine.



A range of endoscopes is available to match different pathologies



Spinal canal decompression with interlaminar access



Intraoperative site after decompression

Today, the instrument sets available permit a full-endoscopic approach under visualization, depending on the indication criteria, which is equivalent to conventional operations. While lateral, stenosis with symptoms on one side can be frequently operated using the basic instrument set, the larger Stenosis System can be used to operate on advanced cases or central stenosis. It is always important to consider whether a stabilizing measure is necessary in addition to decompression.



Priv.-Doz. Dr. med. habil. Sebastian Ruetten



Dr. med. Martin Komp

Center for Spine Surgery and Pain Therapy
Head: Priv.-Doz. Dr. med. habil. Sebastian Ruetten



ST. ELISABETH GRUPPE 
KATHOLISCHE KLINIKEN RHEIN-RUHR

Center for Orthopedics and Traumatology
of the St. Elisabeth Group – Catholic Hospitals Rhein-Ruhr
St. Anna Hospital Herne/Marienhospital Herne University Hospital/Marien Hospital Witten
Director: Prof. Dr. med. Georgios Godolias

VERTEBRIS stenosis


Full-endoscopic, interlaminar decompression

Positioning

The patient is placed on an operating table with an X-ray transparent top in the prone position with pelvic and thorax support pillows. A C-arm image intensifier is required during the procedure.



Prone position with pelvic and thorax pillows



Determination of interlaminar access


Using image intensifier control, the access is determined on the basis of anatomical landmarks in the posterior-anterior beam path and taking account of the pathology. The port must be maximally medial in the interlaminar window in order to permit easier lateral access below the obliquely positioned zygoapophyseal joints.




Marking the entry point on the skin



Entry point should be in a maximally medial position




Access below the zygoapophyseal joints should be possible




Skin incision

Performance of interlaminar access


After determining the entry point on the skin and performance of the skin incision, the dilator is inserted up to the ligamentum flavum or to the zygo-apophyseal joints under posterior-anterior image intensifier control. The subsequent procedure is then performed in the lateral beam path. The working sleeve with oblique opening is pushed over the dilator toward the ligament and the dilator is removed. The endoscope is introduced and the ongoing intervention carried out under continuous visualization and irrigation.



Insertion of the dilator, ...



... operating sleeve ...



... and endoscope

VERTEBRIS stenosis

Full-endoscopic, interlaminar decompression

Ipsilateral, decompression on one side

After the access has been created, the bony structures are exposed. It may be helpful to start decompression at the caudal end of the descending facet. Depending on the pathology, decompression is then commenced with resection of parts of the medial descending facet, the cranial and caudal lamina, and the ligamentum flavum. The extent of decompression generally continues cranially at least until the tip of the ascending facet and caudally to half of the pedicle. The medial portions of the ascending facet and the ligamentum flavum are then resected until sufficient decompression of the neural structures can be clearly seen



It may be helpful to start decompression at the caudal end of the descending facet



The extent of bone resection must generally reach from the tip of the ascending facet to the middle of the caudal pedicle



A range of burrs and bone punches is available for bone resection. They can be introduced through the intraendoscopic working channel



Resection of medial portions of the ascending facet




Removal of protruding annulus and osteophytes



cranially, caudally and laterally. In the case of a central stenosis, the ligamentum flavum generally needs to be resected medially to the midline. Finally, it may be necessary to remove protruding annulus parts and osteophytes in the ventral epidural space. If the patient experiences bilateral symptoms of a lateral stenosis, "over the top" access using the undercutting technique to the opposite side is not carried out. An independent contralateral access is used to retain the median portions of the ligamentum flavum and leave the spinal canal untouched here.



Site after ipsilateral decompression




VERTEBRIS stenosis

Full-endoscopic, interlaminar decompression

Contralateral decompression in over-the-top technique

If bilateral symptoms occur with a central stenosis, a unilateral approach is carried out with "over-the-top" access using the undercutting technique to the opposite side. For this purpose, bone in the ventral area of the spinous process is resected until the contralateral side can be accessed dorsally up to the dura of the spinal cord. If possible, the ligamentum flavum is initially left in place to protect the dura and bony decompression is again carried out by laminotomy and partial facetectomy. The ligamentum flavum is then completely resected. Finally, the contralateral recess needs to be extended. The decompression is completed when the dura and the spinal nerves have been clearly decompressed.




Entrance on the contralateral side

Decompression contralaterally including recess

In general the sealing caps for endoscope and working sleeve should only be used briefly if bleeding obscures visibility since when operations last a long time and the drainage of fluid is prevented without being noticed, the consequences of volume overload and elevated pressure within the spinal canal and the associated and neighbouring structures should not be ignored. An extended and uninterrupted excessive retraction of the neural structures with the working sleeve in a medial direction must be avoided particularly in cranial areas, or only carried out intermittently, in order to avoid the risk of neurological damage. Experience indicates that as with all new techniques there is generally an enhanced risk of problems occurring during the learning curve.



Site after over-the-top decompression





















VERTEBRIS stenosis

Endoscope and accessories

Endoscope	
	<p>Discoscope PANOVIEW Plus Discoscope 20°, working channel ID 5.6 mm, OD 9.3 x 7.4 mm, WL 177 mm (892109205), set incl. sealing cap attachment (8792.452), sealing caps (89.02), sealing membrane (15479.006), O-ring (9500.113) and cleaning brush (6.03) 8921092051</p>
	<p>Endoscope adapter for distance control 892009000</p>
	<p>Fiber light cable D 3.5 mm, WL 1.8m, set including adapter on the projector side (8095.07) and adapter on the endoscope side (809509), color code orange 806635181</p>
	<p>Reprocessing basket for machine reprocessing and sterilization, for discoscopes 89210.xxxx, width 132 mm x length 472 mm x height 74 mm 38044.411</p>
Retainer arm systems	
	<p>Retainer arm adapter for fixing the endoscope in combination with universal retainer arm (898004717) or LEYHLA articulated arm (8766951) 892009070</p>
	<p>Universal retainer arm max. retaining force 90N, joint radius 420 mm 898004717</p>
	<p>Clamp socket electrically insulated to 4.5 KV AC, for attaching to the operating table (standard rail) 8840.9722</p>
	<p>LEYHLA articulated arm for attaching to the operating table, fixing in place using 2 articulated arms is recommended 8766.951</p>










Access and working instruments

Access Instruments	
	Dilator OD 9.4 mm, cannulated, for working sleeve OD 10.5 mm 892209510
Step-dilator Set comprising: 8922095000	
	Dilator OD 3.9 mm, cannulated, for working sleeve 892209505
	Dilator OD 5.9 mm, for working sleeve or dilator OD 7.0 mm 892209507
	Dilator OD 6.9 mm, for working sleeve or dilator OD 9.5 mm 892209508
	Dilator OD 9.4 mm, for working sleeve OD 10.5 mm 892209515
	Working sleeve OD 10.5 mm , WL 120 mm 892209010
	Irrigation adapter OD 10.5 mm 892209310
Working Instruments	
	Tube sheath punches
	Punch , 2.0 mm cutting width, OD 5.5 mm, WL 340 mm, TL 490 mm 892409020
	Punch , 3.5 mm cutting width, OD 5.5 mm, WL 340 mm, TL 490 mm 892409035
	Kerrison punches
	Kerrison punch 60° 4.5 mm cutting width, D 5.5 mm, WL 350 mm, TL 460 mm 892409445
	Kerrison punch 90° 4.5 mm cutting width, D 5.5 mm, WL 350 mm, TL 460 mm 892409945
	Micro punches and rongeurs Color coding for easy identification of instrument diameter
	Rongeur OD 3.0 mm, WL 290 mm 89240.3003
	Rongeur OD 4.0 mm, WL 290 mm 89240.3004
	Punch OD 3.0 mm, WL 290 mm 89240.3023
	Punch OD 4.0 mm, WL 290 mm 89240.3024


VERTEBRIS stenosis

PowerDrive ART1 – Universal motor system








Universal motor system	
Burs for Power Stick M 5	
	Oval burr, with side guard, OD 5.5 mm, WL 290 mm899751505
	Oval burr, eccentric, with side guard, OD 5.5 mm, WL 290 mm899751555
	Round burr, OD 5.5 mm, WL 290 mm899751305
	Round burr, diamond, OD 5.5 mm, WL 290 mm899751405
Articulated burr – TipControl	
	TipControl® – Articulating bone burr, complete, OD 4.0 mm, WL 290 mm (899753754), for Power Stick M5 (15336058), set including cardan burr insert, sterile, pack of 5 (499751704), key for inserting and removing the cardan burr insert (15372005), irrigation adapter (15261106)899753794
Motor handles – Power Stick M5	
	Power Stick M5/0 Handle for shaver blades or burs, operation with footswitch, sterilizable, max. speed 16,000 rpm, with fixed connection cable8995500001
	Power Stick M5/3 Handle for shaver blades or burs, operation with keypad or footswitch, sterilizable, max. speed 16000 rpm, with fixed connection cable8995500031
PowerDrive ART1	
Universal motor system, set incl. power cable, Can Bus connection cable Technical Features: autom. handle and tool recognition, storage function, user-specific parameters, memory function for tools	
	Power supply unit 230 V, 50/60 Hz23040011
	Power supply unit 100 V, 50/60 Hz23040021
	Power supply unit 110 V, 50/60 Hz23040041
	Power supply unit 115 V, 50/60 Hz23040061
	Power supply unit USA 120 V, 50/60 Hz23040071
	Power supply unit 127 V, 50/60 Hz23040121
Power supply unit 240 V, 50/60 Hz23040141	
	Double-pedal footswitch for PowerDrive ART1 (Series 2304)2304.901

COMBIDRIVE EN – High-speed motor system












High-speed motor system

Burrs with distal protection	
	Round burr, tungsten carbide, burr \varnothing 3.0 mm, WL 350 mm, pack of 3 82960.3730
	Support sleeve, with distal protection, OD 4.0 mm 82970.1330
	Round burr, diamond, burr \varnothing 3.0 mm, WL 350 mm, pack of 3 82960.3930
	Support sleeve, with distal protection, AD 4,0 mm 82970.1330
Burrs without distal protection	
	Round burr, diamond, burr \varnothing 3.7 mm, WL 350 mm, pack of 3 82960.3940
	Support sleeve, OD 4.0 mm 82970.1340
High-speed handpiece	
	Handpiece, angled, with adapter 40,000 rpm, INTRA-interface 82950.1301
High-speed motor system COMBIDRIVE EN	
For use with high-speed accessories and accessories for Power Stick M5 (see page 14)	
	COMBIDRIVE EN Set High-speed motor system incl. power cable, footswitch, electronic motor, medium, connection cable and a wide range of accessories 20951.0000
















VERTEBRIS stenosis

Surgitron – Radio frequency unit



RF unit – Surgitron	
	Bipolar accessories – ablation electrodes
	Bipolar hollow sphere electrode distal head part D 2.9 mm, WL 330 mm, disposable899364300
	Bipolar hollow sphere electrode distal head part D 3.4 mm, WL 330 mm, disposable899364400
	Electrode handle, bipolar for mounting bipolar hollow sphere electrodes, fixed cable with US 2 PIN connector, cable length 3 m, reusable899364100
Bipolar accessories – Trigger-Flex	
	Trigger-Flex handle for bipolar controllable Trigger-Flex electrode, including bipolar cable with US 2 PIN connector (2343837)8792.691
	Trigger-Flex bipolar electrode articulated at the distal end, D 2.0 mm, disposable, sterile, pack of 64792.6912
	Trigger-Flex outer sheath WL 400 mm8792.694
	Bipolar connection cable for Trigger-Flex handle, US 2 PIN connector2343.837
RF unit Surgitron	
	Radio frequency unit R.Wolf/elliquence Surgitron Dual 4 MHz, radio frequency unit cpl. with footswitch (2343.901) power for monopolar max. 120 watt/4 MHz, bipolar max. 120 watt/1.7 MHz, 100/120/220/240 V, 50/60 Hz2343.001

Consumables and replacement parts

for TipControl®:	
	Cardanic burr insert, sterile (pack of 5) 499751704
	Wrench 15372005
	Irrigation adapter, complete (M5) 15261106
for Trigger-Flex:	
	Flexible bipolar electrode (pack of 6) 4792.6912
	Trigger-Flex shaft 40 cm 8792.694
	Trigger-Flex replacement seals (pack of 2) 8792.695
	Sealing cap attachment incl. 10 sealing caps (89.00) 8792.452
	Sealing caps opening 0.75 mm for instruments up to D 2.4 mm, black, pack of 10..... 89.00
	Sealing caps opening 2.7 mm for instruments over 3.4 to 5.1 mm, blue, pack of 10..... 89.02
	Sealing membrane 15479006
	Irrigation toggle complete, rotatable..... 15461.034
	O-rings for irrigation toggle (15461.034) pack of 10 9500.113
	Fog reduction agent for endoscopes, disposable, pack of 10 102.02
	Cleaning brush D 5.0 mm, brush length 50 mm, TL 375 mm..... 6.03
	O-rings for irrigation adapter 892209310 15364285

Literature

Ruetten S, Komp M, Hahn P, Oezdemir S.

Decompression of lumbar lateral spinal stenosis: full-endoscopic, interlaminar technique. *Oper Orthop Traumatol* 2013;DOI 10.1007/s00064-012-0195-2

Ruetten S, Komp M, Hahn P, Oezdemir S.

In: Haertl R, Korge A (eds) *Minimally Invasive Spine Surgery – Techniques, Evidence, and Controversies*. Thieme, Stuttgart 2012, pp 59-62

Ruetten S.

Equipment for full-endoscopic spinal surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 59-62

Ruetten S.

Endoscopic lumbar disc surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 303-308

Ruetten S.

Full-endoscopic operations of the spine in disk herniations and spinal stenosis. *Surg Technol Int* 2011;XXI:284-298

Ruetten S.

Full-endoscopic interlaminar lumbar discectomy and spinal decompression. In: Kim DH, Kim K-H, Kim Y-C (eds) *Minimally Invasive Percutaneous Spinal Techniques*. Elsevier, Philadelphia, 2011, pp 351-359

Komp M, Hahn P, Merk H, Godolias G, Ruetten S.

Bilateral operation of lumbar degenerative central spinal stenosis in full-endoscopic interlaminar technique with unilateral approach: Prospective 2-year results of 74 patients. *J Spinal Disord Tech* 2011;24:281-287

Ruetten S, Komp M, Merk H, Godolias G.

Surgical treatment for lumbar lateral recess stenosis with the full-endoscopic interlaminar approach versus conventional microsurgical technique: A prospective, randomized, controlled study. *J Neurosurg Spine* 2009;10:476-485

Ruetten S, Komp M, Merk H, Godolias G.

Recurrent lumbar disc herniation following conventional discectomy: A prospective, randomized study comparing full-endoscopic interlaminar and transforaminal versus microsurgical revision. *J Spinal Disord Tech* 2009;22:122-129

Ruetten S.

Vollendoskopische Operationen der Lendenwirbelsäule. In: Jerosch J, Steinleitner W (ed) *Minimal invasive Wirbelsäulen-Intervention*. Deutscher Ärzte Verlag, Köln, 2009, pp 515-528

Ruetten S, Komp M, Merk H, Godolias G.

Full-endoscopic anterior decompression versus conventional anterior decompression and fusion in cervical disc herniations. *Int Orthop* 2008;33:1677-DOI 10.1007/s00264-008-0684-y

Ruetten S, Komp M, Merk H, Godolias G.

Full-endoscopic cervical posterior foraminotomy for the operation of lateral disc herniations using 5.9-mm endoscopes: A prospective, randomized, controlled study. *Spine* 2008;33:940-948

Ruetten S, Komp M, Merk H, Godolias G.

Full-endoscopic interlaminar and transforaminal lumbar discectomy versus conventional microsurgical technique: A prospective, randomized, controlled study. *Spine* 2008;33:931-939

Ruetten S, Komp M, Merk H, Godolias G.

A new full-endoscopic technique for cervical posterior foraminotomy in the treatment of lateral disc herniations using 6.9-mm endoscopes: prospective 2-year results of 87 patients. *Minim Invas Neurosurg* 2007;50:219-226

Ruetten S, Komp M, Merk H, Godolias G.

Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disc herniations via the interlaminar and lateral transforaminal approach. *J Neurosurg Spine* 2007;6:521-530

Ruetten S, Komp M, Godolias G.

A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6 mm endoscopes: Prospective 2-year results of 331 patients. *Minim Invasive Neurosurg* 2006;49:80-87

Ruetten S, Komp M, Godolias G.

An extreme lateral access for the surgery of lumbar disc herniations inside the spinal canal using the full-endoscopic uniportal transforaminal approach. – Technique and prospective results of 463 patients. *Spine* 2005;30:2570-2578

Ruetten S.

The full-endoscopic interlaminar approach for lumbar disc herniations. In: Mayer HM (ed) *Minimally Invasive Spine Surgery*. Springer, Berlin Heidelberg New York, 2005, pp 346-355

